

Speedcem® Plus

Cementation of root posts (metal, ceramic, FRC)

Befestigung von Wurzelstiften (Metall, Keramik, FRC)


Root Post

Metal / FRC / ceramic: Pre-treatment according to manufacturer's instructions.

Wurzelstift


Metall / FRC / Keramik: Vorbehandlung gemäss Herstellerangaben.

Monobond® Plus



Preparation


Präparation



Cementation

Befestigung


Speedcem® Plus



Core build-up

Stumpaufbau

MultiCore® Flow



more information
www.cementation-navigation.com

Manufacturer
Ivoclar Vivadent AG
9494 Schaan/Liechtenstein
www.ivoclar.com

Rx ONLY
For dental use only!

see
Instructions
www.ivoclar.com/elFU

CE 0123

Date information prepared:
2024-02-20 / Rev. 4

ivoclar