

| Desired 3D MASTER shade | Materials required IPS Style® | | | | | |
|-------------------------------|------------------------------------|-----------------------------------|-----------------------|----------------------------|-----------------------------|-----|
| | IPS Style® Ceram Powder Opaquer | Mixing ratio for the dentin layer | | | IPS Style® Ceram Incisal | |
| OM1 | O BL1/BL2 | D BL2 | | | I BL | |
| OM2 | O BL3/BL4 | D BL2 | | | I BL | |
| OM3 | O BL3/BL4 | D BL3 | | | I BL | |
| 1M1 | O BL1/BL2 | D BL4 | | | I BL | |
| 1M2 | O BL4 | 1/2 D A1 | 1/2 D A2 | | I 1 | |
| 2M1 | O D2 | 2/3 D B1 | 1/3 I 2 | | I 1 | |
| 2L1.5 | O B1 | 1/3 D A2 | 1/3 D B1 | 1/3 D B2 | I 1 | |
| 2R1.5 | O B1 | 4/6 D B1 | 1/6 D D2 | 1/6 CT orange pink | I 1 | |
| 2M2 | O A1 | 1/4 D A1 | 3/4 D A2 | | I 1 | |
| 2L2.5 | O B2 | 1/4 D B3 | 3/4 D B2 | | I 2 | |
| 2R2.5 | O B2 | 4/5 D B2 | 1/5 CT orange pink | | I 2 | |
| 2M3 | O A1 | 1/4 D A2 | 3/4 D B3 | | I 2 | |
| 3M1 | O C1 | 6/10 D2 | 3/10 I 3 | 1/10 CT orange-pink | I 2 | |
| 3L1.5 | O D2 | 3/4 D C1 | 1/4 D C4 | | I 3 | |
| 3R1.5 | O D2 | 4/5 D D2 | 1/5 CT orange-pink | | I 3 | |
| 3M2 | O D3 | 1/2 D A3 | 1/2 D D3 | | I 3 | |
| 3L2.5 | O B3 | 1/2 D B3 | 1/2 D D4 | | I 3 | |
| 3R2.5 | O A3.5 | 3/4 D A3.5 | 1/4 D A2 | | I 3 | |
| 3M3 | O B4 | D B4 | | | I 3 | |
| 4M1 | O C3 | 15/21 D C1 | 2/21 CT khaki | 3/21 SI grey | 1/21 Basic Gingiva BG34 | I 4 |
| 4L1.5 | O C2 | 10/21 D C2 | 10/21 D C3 | 1/21 Basic Gingiva BG34 | | I 4 |
| 4R1.5 | O C2 | 17/20 D C2 | 1/20 SI grey | 2/20 Basic Gingiva BG34 | | I 4 |
| 4M2 | O D3 | 4/6 D C3 | 1/6 CT khaki | 1/6 CT orange | | I 4 |
| 4L2.5 | O A3.5 | 1/4 D A3.5 | 1/4 D C3 | 2/4 D A4 | | I 4 |
| 4R2.5 | O A3.5 | 8/11 D A4 | 1/11 D A3 | 2/11 CT khaki | | I 4 |
| 4M3 | O A3.5 | 5/10 D A4 | 3/10 D B4 | 2/10 CT orange | | I 4 |
| 5M1 | O C3 | 15/20 D C3 | 3/20 SI grey | 2/20 Basic Gingiva BG34 | | I 5 |
| 5M2 | O C3 | 7/10 D A4 | 3/10 CT khaki | | | I 5 |
| 5M3 | O A4 | 6/13 D A4 | 1/13 SI grey | 5/13 OD orange | 1/13 Basic Gingiva BG34 | I 5 |