


Helioseal® F

Fissure sealants don't give the













opponent a chance



Helioseal® F Fissure sealant with power



Step by step

| | | | | | | | | | | | |
|---|---|---|--|---|---|---|---|--|---|--|--|
|  |  |  |  |  |  |  |  |  |  |  |  |
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
| Tooth to be sealed | Clean the surface | Isolate | Etch | Rinse | Dry and check etch pattern | Apply sealant | Light-cure | Check occlusion | Finish | Protected tooth with sealant | Apply fluoride |