

Materials combination table / Massenkombinationstabelle





Chromascop shades / Chromascop Farben

	Bleach	white / weiss	yellow / gelb	light-brown / hellbraun	grey / grau	dark-brown / dunkelbraun
IPS d.SIGN						
Opaquer	 010-020 030-040	 110 120 130 140	 210 220 230 240	 310 320 330 340	 410 420 430 440	 510 520 530 540
Intensiv Opaquer		 white	 violet	 brown	 incisal	
Margin	 010-020 030-040	 110 120 130 140	 210 220 230 240	 310 320 330 340	 410 420 430 440	 510 520 530 540
Intensiv Margin		 yellow	 orange-pink	 add-on	 opaque	 orange
Deep Dentin	 010-020	 110 120 130 140	 210 220 230 240	 310 320 330 340	 410 420 430 440	 510 520 530 540
Dentin	 010-020 030-040	 110 120 130 140	 210 220 230 240	 310 320 330 340	 410 420 430 440	 510 520 530 540
Incisal	 T-S1 T-S1	 S1 S1 S1 S2	 S2 S2 S3 S3	 S3 S3 S3 S3	 S3 S3 S3 S3	 S3 S3 S3 S3
Transparent	 neutral	 neutral	 neutral	 neutral	 neutral	 neutral

Materials combination table / Massenkombinationstabelle

Shades independent of any shade system / Farbsystemunabhängig

IPS d.SIGN Impulse Kit 1				
Occlusal Dentin				
	orange	brown		
Mamelon				
	light	salmon	yellow-orange	
Effect				
	Effect 1	Effect 2	Effect 3	Effect 4
Transparent				
	blue	brown-grey	orange-grey	
Incisal Edge				
	Incisal Edge			
IPS d.SIGN Impulse Kit 2				
Brilliant Dentin				
	orange	white	yellow	
Palatinal Dentin				
	red	yellow		
Insert				
	grey	orange		
Inter Incisal				
	white-blue			
Special Incisal				
	grey	yellow		
Cervical Enamel				
	khaki	orange-pink	orange	yellow

IPS d.SIGN Gingiva Kit								
Gingiva Opaquer								
	Gingiva Opaquer							
Gingiva								
	G1	G2	G3	G4	G5			
Modifier								
	GM1	GM2	GM3	GM4				
IPS d.SIGN Essence Kit								
Essence								
	brown	grey	orange	yellow				
IPS d.SIGN Stains Kit								
Stains								
								
	basic blue	basic red	basic yellow	black	crackliner	grey	khaki	mahogany
	maroon	olive	orange	vanille	white	yellow		
IPS d.SIGN Shade Kit								
Shade								
	1	2	3	4	5	6	7	