

SR Ivocron®

Versatile



PMMA veneering material

One material for three techniques

SR Ivocron is a quality PMMA veneering material for crowns, bridges and temporary restorations. It comes in the 20 shades of the Chromascop shade guide. The material is ideally used with the following Ivoclar Vivadent products: IPS d.SIGN, IPS Empress, SR Chromasit, SR Spectrasit, SR Antaris/ SR Postaris.

The separately sold cervical materials allow the material to be combined with SR Vivodent, SR Vivodent PE, SR Orthotyp and SR Orthotyp PE. SR Ivocron has been successfully used in clinical situations for many years. All in all, it is a very versatile product.



Indication

- Combined denture prosthetics
- Fixed crown and bridge prosthetics
- Temporary technique
- Long-term temporaries for implant prosthetics

Advantages

- Long-term clinical experience
- Three different polymerization techniques
- Small compact assortment
- Universal indications

Pouring

SR Ivocron is mixed with the Cold Liquid for the pouring technique. This technique is predominantly used for temporaries, attaching ground denture teeth to the cast framework and for repairing PMMA.

Layering

SR Ivocron is mixed with the Hot Liquid to achieve a smooth, soft modeling consistency. This material can be used to build up veneers. The material is heat cured under pressure to further enhance its quality.

Pressing

The material is mixed with the Press Liquid to achieve a smooth pressable dough. It is suitable for the flasking technique and can be heat-cured.



Descriptions and data constitute no warranty of attributes and are not binding. The right to make technical changes is reserved. Printed in Austria
© Ivoclar Vivadent AG
563783/0202/e/2/VVA

For more information about the product, please contact Ivoclar Vivadent directly.

Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Fürstentum Liechtenstein
Tel. +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent
technical