

# Programat® S1

## Software Revision Log

Version	Changes from previous version	Status
<b>V6.0</b> 11.03.2019	<p><b>Software optimized</b></p> <ul style="list-style-type: none"> <li>• New program for IPS e.max ZirCAD Prime</li> <li>• Program order modified               <ul style="list-style-type: none"> <li>○ <b>PLEASE NOTE: Individual programs on program number 23 and 24 will be deleted in case of an update!</b></li> </ul> </li> <li>• Various software improvements</li> </ul>	<b>Currently</b> delivered version
<b>V5.0</b> 29.05.2017	<p><b>Software optimized – new program structure</b></p> <ul style="list-style-type: none"> <li>• New programs for IPS e.max ZirCAD MT Multi/MT/LT</li> <li>• Extended function with Programat Dosto-Tray</li> <li>• 50 programs available (22 pre-programmed)               <ul style="list-style-type: none"> <li>○ <b>PLEASE NOTE: All individual programs will be deleted in case of an update!</b></li> </ul> </li> <li>• Various software improvements</li> </ul>	<b>No</b> longer available!
<b>V3.0</b> 09.10.2014	<p><b>Software optimized – new program structure</b></p> <ul style="list-style-type: none"> <li>• New programs for Wieland Dental materials (Zenotec, Zenostar)</li> <li>• Various software improvements</li> </ul>	<b>No</b> longer available!
<b>V2.0</b> 11.10.2012	<p><b>Software optimized – new program structure</b></p> <ul style="list-style-type: none"> <li>• Sintering of IPS e.max ZirCAD single-tooth frameworks in 75 minutes (in the past, 90min were required for this procedure)</li> <li>• Sintering of IPS e.max ZirCAD bridge frameworks with up to 4 units in 90 minutes (in the past, 2h 45min were required for this procedure)</li> <li>• Quick sintering of full-contour ZrO2 single-tooth restorations with enhanced translucency in just 2h 55min</li> <li>• Quick sintering of ZrO2 bridges with enhanced translucency in just 4h 30min</li> </ul>	<b>No</b> longer available!
<b>V1.0</b> 07.12.2009	<b>Basic Version of the Programat S1</b>	<b>No</b> longer available!

Version 11.03.2019

[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Ivoclar Vivadent AG

Bendererstr. 2 | 9494 Schaan | Liechtenstein | Tel.: +423 / 235 35 35 | Fax: +423 / 235 33 60