

**Safety Data Sheet**  
according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

**1 Identification**

- **Product identifier**
- **Trade name:** ***Liquid Strip***
- **Application of the substance / the mixture** *Water-soluble glycerine paste*
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
*Ivoclar Vivadent Inc.*  
*175 Pineview Drive, Amherst, N.Y. 14228*  
*USA*  
*Tel. +1 800 533 6825*  
*Fax +1 716 691 2285*
  
- Ivoclar Vivadent Inc.*  
*1-6600 Dixie Road*  
*Mississauga, Ontario*  
*L5T 2Y2*  
*Canada*  
*Phone: +1 905 670 8499*  
*Fax: +1 905 670 3102*
- **Information department:** *Quality Assurance / Regulatory Affairs*
- **Emergency telephone number:**  
*24 Hour Emergency Assistance:*  
*Emergency-Call USA - Infotrac: 1-800-535-5053*  
*Emergency-Call Canada - Canutec: 1-613-996-6666*
  
- General SDS Assistance:*  
*US: 1-800-533-6825*  
*Canada: 1-800-263-8182*

**2 Hazard identification**

- **Classification of the substance or mixture**  
*The product is not classified according to the Globally Harmonized System (GHS).*

---

- **Label elements**
- **GHS label elements** *Void*
- **Hazard pictograms** *Void*
- **Signal word** *Void*
- **Hazard statements** *Void*
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**
- |            |   |                |
|------------|---|----------------|
| HEALTH     | 0 | Health = 0     |
| FIRE       | 0 | Fire = 0       |
| REACTIVITY | 0 | Reactivity = 0 |

## Safety Data Sheet according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

Trade name: **Liquid Strip**

(Contd. of page 1)

### 3 Composition/Information on ingredients

- **Chemical characterization:** Mixtures
- **Description:** Mixture of the substances listed below with nonhazardous additions.

- **Dangerous components:**

CAS: 56-81-5 glycerol	60-100% w/w
-----------------------	-------------

### 4 First aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:**  
Supply fresh air and to be sure call for a doctor.  
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Rinse with water.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about protection against explosions and fires:** No special measures required.

(Contd. on page 3)

CA

## Safety Data Sheet according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

**Trade name: Liquid Strip**

(Contd. of page 2)

- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
  - Protect from exposure to the light.
  - Protect from heat and direct sunlight.
  - Keep receptacle tightly sealed.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/ Personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

**CAS: 56-81-5 glycerol**

EL Long-term value: 10\* 3\*\* mg/m<sup>3</sup>  
\*mist; \*\*mist, respirable

EV Long-term value: 10 mg/m<sup>3</sup>

- **Additional information:** The lists that were valid during the creation were used as basis.

- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
  - Usual hygienic measures for dental practice and dental laboratories.
  - Wash hands before breaks and at the end of work.
  - Avoid contact with the eyes and skin.
  - Keep away from foodstuffs, beverages and feed.
- **Breathing equipment:** Not required.
- **Protection of hands:**
  - The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
  - Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
  - Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
  - The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
  - The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Safety glasses

### 9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**
  - Form:** Pasty

(Contd. on page 4)

-CA

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

Trade name: **Liquid Strip**

(Contd. of page 3)

· <b>Color:</b>	White
· <b>Odor:</b>	Characteristic
· <b>Odor threshold:</b>	Not determined.
· <b>pH-value:</b>	Not determined.
· <b>Change in condition</b>	
· <b>Melting point/Melting range:</b>	Undetermined.
· <b>Boiling point/Boiling range:</b>	Undetermined.
· <b>Flash point:</b>	Not applicable.
· <b>Auto igniting:</b>	Product is not selfigniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Explosion limits:</b>	
· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.
· <b>Vapor pressure:</b>	Not determined.
· <b>Density:</b>	Not determined.
· <b>Relative density</b>	Not determined.
· <b>Vapor density</b>	Not determined.
· <b>Evaporation rate</b>	Not determined.
· <b>Solubility in / Miscibility with</b>	
· <b>Water:</b>	Nearly insoluble.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
· <b>Dynamic:</b>	Not determined.
· <b>Kinematic:</b>	Not determined.
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid**  
Protect from exposure to the light.  
Keep away from heat and direct sunlight.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.

(Contd. on page 5)

CA

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

**Trade name: Liquid Strip**

(Contd. of page 4)

· **Additional toxicological information:**

*When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.*

· **Carcinogenic categories**

· **NTP (National Toxicology Program)**

*None of the ingredients is listed.*

### 12 Ecological information

· **Toxicity**

· **Aquatic toxicity:** *No further relevant information available.*

· **Persistence and degradability** *No further relevant information available.*

· **Behavior in environmental systems:**

· **Bioaccumulative potential** *No further relevant information available.*

· **Mobility in soil** *No further relevant information available.*

· **Additional ecological information:**

· **General notes:**

*Water hazard class 1 (Self-assessment): slightly hazardous for water*

*Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.*

· **Results of PBT and vPvB assessment**

· **PBT:** *Not applicable.*

· **vPvB:** *Not applicable.*

· **Other adverse effects** *No further relevant information available.*

### 13 Disposal considerations

· **Waste treatment methods**

· **Recommendation:**

*Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.*

· **Uncleaned packagings:**

· **Recommendation:** *Disposal must be made according to official regulations.*

### 14 Transport information

· **UN-Number**

· **DOT, TDG, ADN, IMDG, IATA** *Void*

· **UN proper shipping name**

· **DOT, ADN, IMDG, IATA** *Void*

· **ADR/RID/ADN** *Void*

· **Transport hazard class(es)**

· **DOT, TDG, ADN, IMDG, IATA**

· **Class** *Void*

· **Packing group**

· **DOT, TDG, IMDG, IATA** *Void*

· **Environmental hazards:**

· **Marine pollutant:** *No*

(Contd. on page 6)

CA

## Safety Data Sheet

### according to HPR, Schedule 1

Printing date 05/18/2017

Version number 7

Reviewed on 05/16/2017

**Trade name: *Liquid Strip***

(Contd. of page 5)

· <b>Special precautions for user</b>	<i>Not applicable.</i>
· <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b>	<i>Not applicable.</i>
· <b>Transport/Additional information:</b>	<i>Product is not classified as a dangerous good for transport (ADR, IMDG, IATA).</i>
· <b>UN "Model Regulation":</b>	<i>Void</i>

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

*None of the ingredients is listed.*

· **Section 313 (Specific toxic chemical listings):**

*None of the ingredients is listed.*

· **TSCA (Toxic Substances Control Act):**

*All ingredients are listed.*

· **GHS label elements** *Void*

· **Hazard pictograms** *Void*

· **Signal word** *Void*

· **Hazard statements** *Void*

· **National regulations:**

· **Other regulations, limitations and prohibitive regulations**

*The product is a medical device according to the Directive 93/42/EEC.*

*This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.*

· **Chemical safety assessment:** *A Chemical Safety Assessment has not been carried out.*

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

· **Date of preparation / last revision** 05/18/2017 / 6

· **Abbreviations and acronyms:**

*IMDG: International Maritime Code for Dangerous Goods*

*DOT: US Department of Transportation*

*IATA: International Air Transport Association*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*NFPA: National Fire Protection Association (USA)*

*HMIS: Hazardous Materials Identification System (USA)*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

· **\* Data compared to the previous version altered.**