

Media Release

Schaan, Liechtenstein – 2 December 2013

CAD/CAM-fabricated temporaries

With Telio CAD for Zenotec, Wieland Dental, a company of the Ivoclar Vivadent Group, presents cross-linked PMMA discs for the fabrication of temporary crowns and bridges, as well as implant-retained restorations. Temporaries fabricated of this material are biocompatible and plaque resistant. Moreover, they feature excellent polishability and true-to-nature gloss behaviour.

Telio CAD for Zenotec enables the fabrication of temporary anterior and posterior bridges with up to two consecutive pontics and a maximum wear period of twelve months. It is also suitable for therapeutic restorations in case of corrections of temporomandibular joint problems and the occlusal plane. Given the shade stability and the lifelike fluorescence of the material, a durable, highly esthetic appearance can be achieved. The material is available monochrome in six shades (BL3, A1, A2, A3, A3.5 and B1). The different thicknesses (16 mm and 25 mm) of the discs provide high flexibility of use. Therefore, the material is also suitable for high implant work.

About the Telio System

Telio is a product system for the fabrication of temporary restorations. The materials feature coordinated shades and chemical compatibility and thus offer both dentists and dental technicians high comfort and reliability. Telio CAD for Zenotec restorations are preferably cemented using Telio CS Link or Telio CS Cem Implant after polishing. If desired, the restorations may be individualized with the light-curing Telio Lab LC staining and layering materials or relined using Telio CS C&B.

Telio® is a registered trademark of Ivoclar Vivadent AG.

Zenotec® is a registered trademark of WIELAND Dental + Technik GmbH & Co. KG.

For further information, please contact:

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Phone: +423 235 35 35
Fax: +423 235 33 60
E-mail: info@ivoclarvivadent.com
www.ivoclarvivadent.com

Captions of pictures:

(Telio CAD for Zenotec.jpg)

Fig. 1: Four-unit bridge made of Telio CAD for Zenotec (front) and the shade range.

Media Contact:

Lorenzo Rigliaco
Public Relations Manager
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Principality of Liechtenstein
Phone: +423 235 36 98
Fax: +423 235 36 33
E-mail: lorenzo.rigliaco@ivoclarvivadent.com