

Programat® CS2

Software Revision Log

Version	Changes from previous version	Status
V4.2 11.03.2019	Software optimized <ul style="list-style-type: none"> Optimized parameters for IPS Empress CAD with IPS e.max CAD Crystall. Optimized vacuum test program Optimized dehumidification-function More special characters Acoustic feedback of touchscreen Various languages have been improved Various software optimizations 	Currently delivered version
V4.1 01.02.2018	Software optimized <ul style="list-style-type: none"> Push function can be tested Various languages have been improved Various software optimizations 	No longer available!
V4.0 28.04.2017	Software optimized <ul style="list-style-type: none"> New program structure with new program parameters Various software optimizations 	No longer available!
V3.0 08.04.2016	Software optimized <ul style="list-style-type: none"> In-process dehumidification: Function that protects the vacuum system from contamination and its consequential damages Various software optimizations 	No longer available!
V2.1 09.04.2015	Software optimized <ul style="list-style-type: none"> Various software optimizations 	No longer available!
V2.0 01.10.2014	Software optimized <ul style="list-style-type: none"> Various software optimizations 	No longer available!

Programat® CS2

V1.1 06.06.2014	Software optimized <ul style="list-style-type: none"> The software was attuned to meet the requirements of the Programat App (for iPhones) and various software optimizations implemented. The software was optimized to meet the requirements of the PrograBase X10 Software. As a result, more functions with regard to PrograBase X10 are available with the Programat CS2. 	No longer available!
V1.0 20.09.2013	Basic version	No longer available!

Version 11.03.2019