

## Media Release

Schaan, Liechtenstein – 20 October 2016

---

# Zenostar: Now even greater esthetics and efficiency

## The portfolio extension with pre-shaded discs offers new possibilities

The pre-shaded Zenostar MT discs in shades of MT 1 to 4 are suitable for the efficient production of zirconium oxide restorations. Shade infiltration and pre-drying is not required. This saves valuable time.

Zenostar MT 1 - 4 complement the existing product range. Their level of translucency is equivalent to that of the non-shaded Zenostar MT 0. This coordinated light transmittance is characterized by a natural, esthetic appearance in the patient's mouth. The shade concept is designed so that restorations in the most popular A - D shades only require one stain and one glaze bake after sintering. The desired tooth shade can therefore be reproduced more easily and efficiently.

### Two different disc thicknesses

With Zenostar MT 1- 4 it is possible to produce full-contour single crowns and up to three-unit bridges in the anterior and posterior region. The new discs are available in thicknesses of 14 mm and 18 mm. Shade adaptations or individualized characterizations can be realized with the stains from IPS Ivocolor range.

### Reliable and quick cementation

The self-curing resin cement SpeedCEM® Plus is recommended for the cementation of Zenostar MT restorations in the patient's mouth. Due to the cement's self-adhesive properties, the restorations do not need to be conditioned and the application of a primer is unnecessary. This also increases efficiency.

*Zenostar® is a registered trademark of Wieland Dental + Technik GmbH & Co. KG.  
SpeedCEM® is a registered trademark of Ivoclar Vivadent AG.*

Captions:

*(Zenostar.jpg)*

Fig. More efficiency and esthetics: Zenostar MT 1 - 4

**For more information:**

Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
Tel.: +423 / 235 35 35  
Fax: +423 / 235 33 60  
E-Mail: [info@ivoclarvivadent.com](mailto:info@ivoclarvivadent.com)  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Media Contact:**

André Büssers  
Public Relations Manager  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein  
Tel.: +423 / 235 36 98  
Fax: +423 / 235 36 33  
E-Mail: [andre.buessers@ivoclarvivadent.com](mailto:andre.buessers@ivoclarvivadent.com)