

PrograMill PM3/5/7

PrograMill Software Updates V86.89 (PM7) / V6.26 (PM5) / V4.25 (PM3) *Changelog*

New features and enhancements that apply to all machines:

- **new feature:** Update translations for English, French, Italian and Korean
- **new feature:** Material changer automatically turns to next free position
- **new feature:** New default images for all new materials
- **new feature:** Add information of TeamViewer ID in Machine Data tab
- **new feature:** Add information of OS Version in the Settings tab
- **new feature:** Change IoT default settings, IoT is now on per default
- **new feature:** It is now possible to have multiple Non-Ivoclar Digital tools in one magazine
- **new feature:** Repeat a job from the beginning or to continue the job from the current block
- **new feature:** Grey out stop button during disc load / unload process
- **new feature:** Enhanced Tool Line behavior, pick and place tool directly from Tool Line into magazine
- **new feature:** B-Axis behavior when front flap is open, show animation in case B-Axis is out of range
- **new feature:** New filter for Non-Ivoclar Digital tools in the Tool Library
- **new feature:** Show software version and serial number of the machine in error screens and popups
- **new feature:** Show navigation arrows in error screen in case of multiple errors
- **new feature:** New Cleaning Concept with additional cleaning processes
- Several minor and critical bug fixes