

## Media Release

Cologne (D) / Schaan (FL) – 21 March 2011

---

### Efficient restoration of posterior teeth

**Tetric EvoCeram Bulk Fill is the new composite for the direct restoration of posterior teeth. This material enables cavities to be filled with increments of up to 4 mm.**

Tetric EvoCeram Bulk Fill from Ivoclar Vivadent is a mouldable posterior composite suitable for the bulk-filling technique. Due to its special consistency, Tetric EvoCeram Bulk Fill readily adapts to the cavity walls. The composite is based on the same chemistry as the clinically proven, light-curing, universal nano-hybrid composite Tetric EvoCeram, which is used successfully around the world.

#### **One increment, one filling, one material**

Complicated incremental layering is not required with Tetric EvoCeram Bulk Fill. Instead, an efficient and easy placement technique can be employed, as the composite enables bulk-filling of up to 4 mm. The decisive advantage: No additional universal composite is required to create a covering layer in the enamel area. The Turbo mode of the bluephase 20i curing light (2000 mW/cm<sup>2</sup>) allows composite increments of up to 4 mm to be cured in five seconds. Moreover, TetricEvoCeram Bulk Fill can be anatomically shaped and contoured, so that finishing time is minimized. Because of the well-balanced filler composition, high-gloss polishing of Tetric EvoCeram Bulk Fill restorations is quick and easy to accomplish.

#### **Natural shade adaptation**

Tetric EvoCeram Bulk Fill is available in three universal shades. Its harmonious integration with the natural tooth structure provides the basis for the efficient creation of restorations that closely mimic their natural counterpart. The composite is suitable for the quick restoration of deciduous teeth as well as the efficient restoration of the permanent dentition.

Tetric EvoCeram Bulk Fill will be available from May 2011 onwards.

Tetric EvoCeram® and bluephase® are registered trademarks of Ivoclar Vivadent AG.

For further information please contact:

Ivoclar Vivadent AG  
Bendererstr. 2  
9494 Schaan  
Principality of Liechtenstein  
Tel.: +423 235 35 35  
Fax: +423 235 33 60  
E-mail: [info@ivoclarvivadent.com](mailto:info@ivoclarvivadent.com)  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Caption:

*(Tetric EvoCeram Bulk Fill.jpg)*

Fig. 1: The new Tetric EvoCeram Bulk Fill composite from Ivoclar Vivadent.

---

Media contact:

Lorenzo Rigliaco  
Public Relations Manager  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan  
Principality of Liechtenstein  
Tel.: +423 235 36 98  
Fax: +423 235 36 33  
E-mail: [lorenzo.rigliaco@ivoclarvivadent.com](mailto:lorenzo.rigliaco@ivoclarvivadent.com)