

Media Release

Cologne (D) / Schaan (FL) – 21 March 2011

Semi-permanently cementing restorations on implant abutments

Telio CS Cem Implant is the new dual-curing resin cement from Ivoclar Vivadent. It has been designed for the semi-permanent cementation of restorations on implant abutments.

Telio is a comprehensive temporary restoration system which is geared to dentists, CAD/CAM users and dental technicians alike. In the course of the further development of the Telio range, Telio CS Cem Implant has been introduced, a resin cement designed for the semi-permanent cementation of restorations on implant abutments.

Esthetic cementation of long-term temporaries

Telio CS Cem Implant is suitable for cementing metal-borne and metal-free long-term temporaries made of e.g Telio Lab or Telio CAD. In addition, the material can be used for the semi-permanent cementation of permanent crowns and bridges on implant abutments. Due to the translucency of the cement, highly esthetic results are achieved. Telio CS Cem Implant has been designed for wear periods of six months and more. The cement is free of eugenol and it is radiopaque. It is supplied in a double-push syringe which automatically mixes the components and allows the precise amount of cement to be dispensed. Due to static mixing tip, freshly and homogenously mixed material is available for each application.

Telio CS Cem Implant will be available in the shades Transparent and Opaque from April 2011 onwards.

Telio® is a registered trademark of Ivoclar Vivadent AG.

For further information please contact:

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Principality of Liechtenstein
Tel.: +423 235 35 35
Fax: +423 235 33 60
E-mail: info@ivoclarvivadent.com
www.ivoclarvivadent.com

Caption:

(Telio CS Cem Implant.jpg)

Fig. 1: Telio CS Cem Implant: The new dual-curing resin cement from Ivoclar Vivadent

Media contact:

Lorenzo Rigliaco
Public Relations Manager
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Principality of Liechtenstein
Tel.: +423 235 36 98
Fax: +423 235 36 33
E-mail: lorenzo.rigliaco@ivoclarvivadent.com